

GENESIS SCLERAL

“Created from Proven Technology”

M C T H

0.2

A F D Z E

0.3

B G L Y C K I

0.4

H D F Z V X T

0.5

D L V B N C M F

0.6

Custom Contact Lenses




800-223-1858
www.metro-optics.com

GENESIS SCLERAL

Born from the most proven technology in the contact lens industry, Genesis Scleral provides a powerful and dependable system for addressing irregular corneas where other modalities cannot provide a stable fit or good visual acuity.

- Proven Technology
- Versatile Application
- Stable & Dependable
- Expert Technical Support
- Rapid Delivery

Genesis Scleral can be used on a wide range of corneal abnormalities and extreme dry eye. Genesis Scleral is available in two main series: Genesis 156 (mini-scleral) with a 15.6 mm OAD and Genesis 182 with an 18.2 mm OAD. Both of the Genesis Scleral



Topography of Irregular Cornea

lens versions vault over the cornea landing on the scleral region of the eye and are designed to create a significant tear pool around the cornea. Patients wearing Genesis Scleral experience excellent Visual Acuity, Comfort and Corneal Health.

Genesis Scleral is easily prescribed using a fitting kit. The diagnostic kits feature a versatile selection of lenses that allow the design to be fit on a wide range of cases including: keratoconus, pellucid marginal degeneration, corneal transplant, post-refractive surgery, post trauma and severe dry eye. Genesis Scleral is also backed by superior service & support with decades of experience. The Metro Optics customer service team is friendly, knowledgeable and ready to process orders quickly. Our technical support team of highly experienced consultants is always standing by to assist with design and fitting.



Warranty

Genesis Scleral is warranted for 90 days from initial order. It may be remade as often as necessary and/or returned for full refund of purchase price within the warranty period.

Return & Credit Policy

Genesis Scleral lenses returned for refund must be in factory packaging & accompanied by original invoicing. Credit is not issued on shipping & handling fees. Return shipping fees and proof of delivery on all returned product is the responsibility of the customer.