

Custom Contact Lenses



Custom Soft Lenses for the Surgically Altered Cornea!

RevitalEyes[®]

P O S T S U R G I C A L

Easy to Fit! RevitalEyes features a simple diagnostic fitting system for a fast, easy and accurate fit!


RevitalEyes is a **custom manufactured soft lens** specifically designed for the surgically altered cornea where visual acuity is less than optimal. It is available in two lens material options: standard hydrogel Hioxifilcon B 49% water or Definitive SiHy 74% water material.

Note: for Corneal Graft cases, only Definitive SiHy should be used.

RevitalEyes is as easy to use as a standard commodity planned replacement soft lens!

Contact Metro Optics today at 800-223-1858 to learn more and to order your Custom Fitting Kit. Learn more online at www.metro-optics.com.

This lens is FDA cleared for PKP as well as post refractive surgery.



Actual Corneal Topography of a Post Refractive Surgery Patient



Order Your Fitting Kit Today