



## MANAGE MYOPIA THE EASY WAY!

The Ortho-ATX is a custom designed overnight orthokeratology lens that utilizes an empiric fitting approach to treat high amounts of myopia and astigmatism, producing highly successful outcomes quickly and easily.

- **No Fitting Set Required - Custom lenses designed from spectacle Rx and corneal topography measurements.**
- **Free online training**
- **90% first lens success rate**

[metro-optics.com](http://metro-optics.com)


15802 Vision Drive  
Pflugerville, Texas 78660  
(800) 223-1858

[info@metro-optics.com](mailto:info@metro-optics.com)





Ortho-ATX features an optimized four zone configuration comprised of unique proprietary shaping curves and peripheral alignment stabilization curves to keep the Ortho-ATX centered directly over the desired treatment area of the cornea for highly predictable and reliable results.



- **Effective on Myopia up to -5.00 Diopters**
- **Treats significant amounts of corneal astigmatism**
- **Reverse curve shaping and recovery zones work together to provide desired corneal molding.**
- **Alignment Curve Stabilization provides a stable and comfortable alignment fit across the mid-peripheral cornea while providing adequate tear exchange.**



Vision Shaping  
Treatment VST®

[metro-optics.com](http://metro-optics.com)