



# METRO **K**

for Ortho-K

## MANAGE MYOPIA THE EASY WAY!

The Metro-K is a custom designed overnight orthokeratology lens that utilizes an empiric fitting approach to treat high amounts of myopia and astigmatism, producing highly successful outcomes quickly and easily.

- **No Fitting Set Required — Custom contact lenses designed from spectacle Rx and corneal topography measurements**
- **Free online training**
- **90% first lens success rate**

[metro-optics.com](http://metro-optics.com)

15802 Vision Drive  
Pflugerville, Texas 78660  
(800) 223-1858


[info@metro-optics.com](mailto:info@metro-optics.com)

Custom Contact Lenses  
**METRO**  
**OPTICS**

# METRO K

for Ortho-K

Metro-K features an optimized four zone configuration comprised of unique proprietary shaping curves and peripheral alignment stabilization curves to keep the Metro-K centered directly over the desired treatment area of the cornea for highly predictable and reliable results.



- **Effective on Myopia up to -6.00 Diopters**
- **Treats up to 1.75 Diopters of corneal astigmatism**
- **Reverse curve shaping and recovery zones work together to provide desired corneal molding.**
- **Alignment Curve Stabilization provides a stable and comfortable alignment fit across the mid-peripheral cornea while providing adequate tear exchange**



**OPTIMUM**  
**INFINITE**  
ORTHO-K

[metro-optics.com](http://metro-optics.com)